

PACO28570

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of a putative tumor-recognition complex. Involved in the function of NK cells.

Gene ID:

NKTR

Uniprot

P30414

Synonyms:

NK-tumor recognition protein (NK-TR protein) (Natural-killer cells cyclophilin-related protein) [Includes: Putative peptidyl-prolyl cis-trans isomerase (PPIase) (EC 5.2.1.8) (Rotamase)], NKTR

Immunogen:

Recombinant Human NK-tumor recognition protein (10-175AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A