

# NR2C2AP Antibody, HRP conjugated



PACO28643

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May act as a repressor of NR2C2-mediated transactivation by suppressing the binding between NR2C2/TR4 and the TR4-response element in target genes.

**Gene ID:**

NR2C2AP

**Uniprot**

Q86WQ0

**Synonyms:**

Nuclear receptor 2C2-associated protein (TR4 orphan receptor-associated 16 kDa protein), NR2C2AP, TRA16

**Immunogen:**

Recombinant Human Nuclear receptor 2C2-associated protein (7-139AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A