

# NRARP Antibody, FITC conjugated



PACO28648

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Downstream effector of Notch signaling. Involved in the regulation of liver cancer cells self-renewal. May play a role in the formation of somites.

**Gene ID:**

NRARP

**Uniprot**

Q7Z6K4

**Synonyms:**

Notch-regulated ankyrin repeat-containing protein, NRARP

**Immunogen:**

Recombinant Human Notch-regulated ankyrin repeat-containing protein (1-114AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A