

NTS Antibody, FITC conjugated



PACO28656

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Neurotensin may play an endocrine or paracrine role in the regulation of fat metabolism. It causes contraction of smooth muscle.

Gene ID:

NTS

Uniprot

P30990

Synonyms:

Neurotensin/neuromedin N [Cleaved into: Large neuromedin N (NmN-125); Neuromedin N (NN) (NmN); Neurotensin (NT); Tail peptide], NTS

Immunogen:

Recombinant Human Neurotensin/neuromedin N protein (24-143AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A