

# NTS Antibody, Biotin conjugated



PACO28657

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Neurotensin may play an endocrine or paracrine role in the regulation of fat metabolism. It causes contraction of smooth muscle.

**Gene ID:**

NTS

**Uniprot**

P30990

**Synonyms:**

Neurotensin/neuromedin N [Cleaved into: Large neuromedin N (NmN-125); Neuromedin N (NN) (NmN); Neurotensin (NT); Tail peptide], NTS

**Immunogen:**

Recombinant Human Neurotensin/neuromedin N protein (24-143AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A