

OTX2 Antibody, Biotin conjugated



PACO28773

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probably plays a role in the development of the brain and the sense organs. Can bind to the BCD target sequence (BTS): 5'-TCTAATCCC-3'.

Gene ID:

OTX2

Uniprot

P32243

Synonyms:

Homeobox protein OTX2 (Orthodenticle homolog 2), OTX2

Immunogen:

Recombinant Human Homeobox protein OTX2 protein (1-297AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A