

# PAH Antibody, Biotin conjugated



PACO28797

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

cytosol, extracellular exosome, phenylalanine 4-monooxygenase activity, catecholamine biosynthetic process, cellular amino acid, biosynthetic process, cellular nitrogen compound metabolic process, L-phenylalanine catabolic process, neurotransmitter biosynthetic process, small molecule metabolic process.

**Gene ID:**

PAH

**Uniprot**

P00439

**Synonyms:**

Phenylalanine-4-hydroxylase (PAH) (EC 1.14.16.1) (Phe-4-monooxygenase), PAH

**Immunogen:**

Recombinant Human Phenylalanine-4-hydroxylase protein (2-452AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A