

PGM1 Antibody, FITC conjugated



PACO28908

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This enzyme participates in both the breakdown and synthesis of glucose.

Gene ID:

PGM1

Uniprot

P36871

Synonyms:

Phosphoglucomutase-1 (PGM 1) (EC 5.4.2.2) (Glucose phosphomutase 1), PGM1

Immunogen:

Recombinant Human Phosphoglucomutase-1 protein (1-562AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A