

PHF6 Antibody, HRP conjugated



PACO28927

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional regulator that associates with ribosomal RNA promoters and suppresses ribosomal RNA (rRNA) transcription.

Gene ID:

PHF6

Uniprot

Q8IWS0

Synonyms:

PHD finger protein 6 (PHD-like zinc finger protein), PHF6, CENP-31 KIAA1823

Immunogen:

Recombinant Human PHD finger protein 6 protein (2-312AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A