

PACO28938

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

PIGR

Uniprot

P81265

Synonyms:

Polymeric immunoglobulin receptor (PIgR) (Poly-Ig receptor) [Cleaved into: Secretory component], PIGR


Immunogen:

Recombinant Bovine Polymeric immunoglobulin receptor protein (399-599AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Polymeric immunoglobulin receptor antibody at 2 μ g/ml + HL60 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 83, 59 kDa

Observed band size: 83 kDa

.