

PACO28986

Product Information

Size:

50ug

Reactivity:

Rat, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Acts as a receptor for urokinase plasminogen activator. Plays a role in localizing and promoting plasmin formation. Mediates the proteolysis-independent signal transduction activation effects of U-PA.

Gene ID:

Plaur

Uniprot

P49616

Synonyms:

Urokinase plasminogen activator surface receptor (U-PAR) (uPAR) (CD_antigen: CD87), Plaur


Immunogen:

Recombinant Rat Urokinase plasminogen activator surface receptor ptotein (25-299AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Plaur antibody at 2 μ g/ml

Lane 1: 293T whole cell lysate

Lane 2: Rat brain tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 36, 25 kDa

Observed band size: 36 kDa

.