

# PLP1 Antibody, FITC conjugated



PACO29014

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

This is the major myelin protein from the central nervous system. It plays an important role in the formation or maintenance of the multilamellar structure of myelin.

**Gene ID:**

PLP1

**Uniprot**

P60201

**Synonyms:**

Myelin proteolipid protein (PLP) (Lipophilin), PLP1, PLP

**Immunogen:**

Recombinant Human Myelin proteolipid protein (143-197AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A