

PLP1 Antibody, Biotin conjugated



PACO29016

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This is the major myelin protein from the central nervous system. It plays an important role in the formation or maintenance of the multilamellar structure of myelin.

Gene ID:

PLP1

Uniprot

P60201

Synonyms:

Myelin proteolipid protein (PLP) (Lipophilin), PLP1, PLP

Immunogen:

Recombinant Human Myelin proteolipid protein (143-197AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A