

PMP22 Antibody



PACO29026

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Might be involved in growth regulation, and in myelination in the peripheral nervous system.

Gene ID:

PMP22

Uniprot

Q01453

Synonyms:

Peripheral myelin protein 22 (PMP-22) (Growth arrest-specific protein 3) (GAS-3), PMP22, GAS3


Immunogen:

Recombinant Human Peripheral myelin protein 22 protein (32-64AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO29026 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).