

PRAMEF5 Antibody



PACO29152

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

PRAMEF5

Uniprot

Q5TYX0

Synonyms:

PRAME family member 5, PRAMEF5, PRAMEF23 PRAMEF5L

Immunogen:

Recombinant Human PRAME family member 5 protein (1-476AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A