

Trypsin Antibody



PACO29172

Product Information

Size:

50ug

Reactivity:

Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P00761

Synonyms:

Trypsin (EC 3.4.21.4)

Immunogen:

Recombinant Pig Trypsin protein (9-231AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A