

PACO29208

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular space, serine-type endopeptidase activity, digestion, proteolysis, response to caffeine, response to nicotine, response to nutrient, response to organic substance.

Gene ID:

Prss1

Uniprot

P00762

Synonyms:

Anionic trypsin-1 (EC 3.4.21.4) (Anionic trypsin I) (Pretrypsinogen I) (Serine protease 1), Prss1, Try1

Immunogen:

Recombinant Rat Anionic trypsin-1 protein (24-246AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A