

PACO29244

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the dephosphorylation of phosphoserine (P-Ser). Also catalyzes the hydrolysis of phosphothreonine (P-Thr).

Gene ID:

serB

Uniprot

P0AGB0

Synonyms:

Phosphoserine phosphatase (PSP) (PSPase) (EC 3.1.3.3) (O-phosphoserine phosphohydrolase), serB

Immunogen:

Recombinant Escherichia coli Phosphoserine phosphatase protein (1-322AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A