

RAB2A Antibody, FITC conjugated



PACO29326

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for protein transport from the endoplasmic reticulum to the Golgi complex.

Gene ID:

RAB2A

Uniprot

P61019

Synonyms:

Ras-related protein Rab-2A, RAB2A, RAB2

Immunogen:

Recombinant Human Ras-related protein Rab-2A protein (2-212AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A