

RBP5 Antibody



PACO29416

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Intracellular transport of retinol.

Gene ID:

RBP5

Uniprot

P82980

Synonyms:

Retinol-binding protein 5 (Cellular retinol-binding protein III) (CRBP-III) (HRBPiso), RBP5


Immunogen:

Recombinant Human Retinol-binding protein 5 protein (1-135AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: RBP5 antibody at 8 μ g/ml + Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 16 kDa

Observed band size: 20 kDa

.