

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Bactericidal C-type lectin. Regulates keratinocyte proliferation and differentiation after skin injury via activation of EXTL3-PI3K-AKT signaling pathway.

Gene ID:

Reg3a

Uniprot

P35231

Synonyms:

Regenerating islet-derived protein 3-alpha (REG-3-alpha) (Islet of Langerhans regenerating protein 3) (REG 3) (Lithostathine 3) (Pancreatitis-associated protein 2) (RegIII) (Regenerating islet-derived protein III-alpha) (Reg III-alpha) [Cleaved into: Regenerating islet-derived protein 3-alpha 16.5 kDa form; Regenerating islet-derived protein 3-alpha 15 kDa form], Reg3a, Pap2 Reg3

Immunogen:

Recombinant Rat Regenerating islet-derived protein 3- α ; protein (26-174AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A