

RFPL4B Antibody



PACO29472

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

RFPL4B

Uniprot

Q6ZWI9

Synonyms:

Ret finger protein-like 4B (RING finger protein 211), RFPL4B, RNF211


Immunogen:

Recombinant Human Ret finger protein-like 4B protein (1-263AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO29472 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO29472 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human colon cancer using PACO29472 at dilution of 1:100.