

RNF141 Antibody



PACO29536

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in spermatogenesis.

Gene ID:

RNF141

Uniprot

Q8WVD5

Synonyms:

RING finger protein 141 (Zinc finger protein 230), RNF141, ZNF230

Immunogen:

Recombinant Human RING finger protein 141 protein (1-230AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A