

ZNF91 Antibody



PACO29576

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

ZNF91

Uniprot

Q05481

Synonyms:

Zinc finger protein 91 (Zinc finger protein HPF7) (Zinc finger protein HTF10), ZNF91


Immunogen:

Recombinant Human Zinc finger protein 91 protein (1-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO29576 at dilution of 1:100.