

ROPN1B Antibody



PACO29584

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

ROPN1B

Uniprot

Q9BZX4

Synonyms:

Ropporin-1B (Rhopilin-associated protein 1B), ROPN1B


Immunogen:

Recombinant Human Ropporin-1B protein (1-120AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO29584 at dilution of 1:100.