

RPL34 Antibody



PACO29604

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm, cytosol, cytosolic large ribosomal subunit, extracellular exosome, nucleolus, RNA binding, structural constituent of ribosome, cellular protein metabolic process, gene expression, nuclear-transcribed mRNA catabolic process, nonsense-mediated decay.

Gene ID:

RPL34

Uniprot

P49207

Synonyms:

60S ribosomal protein L34 (Large ribosomal subunit protein eL34), RPL34


Immunogen:

Recombinant Human 60S ribosomal protein L34 protein (2-117AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO29604 at dilution of 1:100.