

S100A11 Antibody



PACO29684

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Facilitates the differentiation and the cornification of keratinocytes.

Gene ID:

S100A11

Uniprot

P31949

Synonyms:

Protein S100-A11 (Calgizzarin) (Metastatic lymph node gene 70 protein) (MLN 70) (Protein S100-C) (S100 calcium-binding protein A11) [Cleaved into: Protein S100-A11, N-terminally processed], S100A11, MLN70 S100C


Immunogen:

Recombinant Human Protein S100-A11 protein (2-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO29684 at dilution of 1:100.