

SAA1 Antibody



PACO29736

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Major acute phase reactant. Apolipoprotein of the HDL complex.

Gene ID:

SAA1

Uniprot

P19708

Synonyms:

Serum amyloid A protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA)], SAA1

Immunogen:

Recombinant Dog Serum amyloid A protein (19-129AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A