

SCGB2A1 Antibody



PACO29752

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May bind androgens and other steroids, may also bind estramustine, a chemotherapeutic agent used for prostate cancer. May be under transcriptional regulation of steroid hormones.

Gene ID:

SCGB2A1

Uniprot

O75556

Synonyms:

Mammaglobin-B (Lacryglobin) (Lipophilin-C) (Mammaglobin-2) (Secretoglobin family 2A member 1), SCGB2A1, LIPHC MGB2 UGB3

Immunogen:

Recombinant Human Mammaglobin-B protein (19-95AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A