

Serpina1 Antibody



PACO29828

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibitor of serine proteases. The primary target is elastase, but also has a moderate affinity for plasmin and thrombin.

Gene ID:

Serpina1

Uniprot

P17475

Synonyms:

Alpha-1-antiproteinase (Alpha-1-antitrypsin) (Alpha-1-proteinase inhibitor) (Serpina A1), Serpina1

Immunogen:

Recombinant Rat Alpha-1-antiproteinase protein (25-411AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A