

Serpib2 Antibody, HRP conjugated



PACO29841

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Inhibits urokinase-type plasminogen activator.

Gene ID:

Serpib2

Uniprot

P29524

Synonyms:

Plasminogen activator inhibitor 2 type A (PAI-2) (PAI2A) (Serpib B2), Serpib2, Pai2, Planh2

Immunogen:

Recombinant Rat Plasminogen activator inhibitor 2 type A protein (1-416AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A