

# Serpinf1 Antibody



PACO29848

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Serpinf1

**Uniprot**

Q80ZA3

**Synonyms:**

Pigment epithelium-derived factor (PEDF) (Caspin) (Serpin F1) (Stromal cell-derived factor 3) (SDF-3), Serpinf1, Pedf Sdf3

**Immunogen:**

Recombinant Rat Alpha-2 antiplasmin protein (20-418AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A