

PACO29868

---

## Product Information

**Size:**

50ug

**Reactivity:**

Bovine

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Pulmonary surfactant associated proteins promote alveolar stability by lowering the surface tension at the air-Liquid interface in the peripheral air spaces.

**Gene ID:**

SFTPC

**Uniprot**

P15783

**Synonyms:**

Pulmonary surfactant-associated protein C (SP-C) (Pulmonary surfactant-associated proteolipid SPL(Val)), SFTPC, SFTP2

**Immunogen:**

Recombinant Bovine Pulmonary surfactant-associated protein C protein (25-58AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A