

SFTPC Antibody, HRP conjugated



PACO29871

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Pulmonary surfactant associated proteins promote alveolar stability by lowering the surface tension at the air-Liquid interface in the peripheral air spaces.

Gene ID:

SFTPC

Uniprot

P11686

Synonyms:

Pulmonary surfactant-associated protein C (SP-C) (Pulmonary surfactant-associated proteolipid SPL(Val)) (SP5), SFTPC, SFTP2

Immunogen:

Recombinant Human Pulmonary surfactant-associated protein C protein (24-58AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A