

# SLITRK3 Antibody



PACO29984

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Suppresses neurite outgrowth.

**Gene ID:**

SLITRK3

**Uniprot**

O94933

**Synonyms:**

SLIT and NTRK-like protein 3, SLITRK3, KIAA0848

**Immunogen:**

Recombinant Human SLIT and NTRK-like protein 3 protein (746-977AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A