

Snca Antibody



PACO30020

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in the regulation of dopamine release and transport.

Gene ID:

Snca

Uniprot

P37377

Synonyms:

Alpha-synuclein, Snca

Immunogen:

Recombinant Rat Alpha-synuclein protein (1-140AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A