

Sod1 Antibody



PACO30072

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Destroys radicals which are normally produced within the cells and which are toxic to biological systems.

Gene ID:

Sod1

Uniprot

P08228

Synonyms:

Superoxide dismutase [Cu-Zn] (EC 1.15.1.1), Sod1

Immunogen:

Recombinant Mouse Superoxide dismutase [Cu-Zn] protein (2-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A