

# SOD2 Antibody, Biotin conjugated



PACO30079

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Destroys superoxide anion radicals which are normally produced within the cells and which are toxic to biological systems.

**Gene ID:**

SOD2

**Uniprot**

P04179

**Synonyms:**

Superoxide dismutase [Mn], mitochondrial (EC 1.15.1.1), SOD2

**Immunogen:**

Recombinant Human Superoxide dismutase [Mn], mitochondrial protein (64-175AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A