

# SRP54 Antibody, Biotin conjugated



PACO30135

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds to the signal sequence of presecretory protein when they emerge from the ribosomes and transfers them to TRAM (translocating chain-associating membrane protein).

**Gene ID:**

SRP54

**Uniprot**

P61011

**Synonyms:**

Signal recognition particle 54 kDa protein (SRP54), SRP54

**Immunogen:**

Recombinant Human Signal recognition particle 54 kDa protein (1-504AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A