

Sry Antibody



PACO30136

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Sry

Uniprot

Q05738

Synonyms:

Sex-determining region Y protein (Testis-determining factor), Sry, Tdf Tdy

Immunogen:

Recombinant Mouse Sex-determining region Y protein (1-143AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A