

PACO30242

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Functions in exocytosis in pancreatic acinar cells regulating the fusion of zymogen granules with each other. May have a pore-forming activity on membranes and regulate exocytosis in other exocrine tissues.

Gene ID:

SYCN

Uniprot

Q0VAF6

Synonyms:

Syncollin (Insulin synthesis-associated protein 1), SYCN, INSSA1 SYL


Immunogen:

Recombinant Human Syncollin protein (22-134AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO30242 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO30242 at dilution of 1:100.