

TARS2 Antibody, HRP conjugated



PACO30267

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

gene expression, mitochondrial threonyl-tRNA aminoacylation, tRNA aminoacylation for protein translation.

Gene ID:

TARS2

Uniprot

Q9BW92

Synonyms:

Threonine--tRNA ligase, mitochondrial (EC 6.1.1.3) (Threonyl-tRNA synthetase) (ThrRS) (Threonyl-tRNA synthetase-like 1), TARS2, TARSL1

Immunogen:

Recombinant Human Threonine--tRNA ligase, mitochondrial protein (369-718AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A