

TBX18 Antibody



PACO30298

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transcriptional regulator involved in developmental processes.

Gene ID:

TBX18

Uniprot

O95935

Synonyms:

T-box transcription factor TBX18 (T-box protein 18), TBX18

Immunogen:

Recombinant Human T-box transcription factor TBX18 protein (1-607AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A