

Thbs1 Antibody



PACO30370

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. Binds heparin. May play a role in dentinogenesis and/or maintenance of dentin and dental pulp. Ligand for CD36 mediating antiangiogenic properties.

Gene ID:

Thbs1

Uniprot

P35441

Synonyms:

Thrombospondin-1, Thbs1, Tsp1

Immunogen:

Recombinant Mouse Thrombospondin-1 protein (19-350AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A