

Thbs3 Antibody



PACO30378

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. Can bind to fibrinogen, fibronectin, laminin and type V collagen.

Gene ID:

Thbs3

Uniprot

Q05895

Synonyms:

Thrombospondin-3, Thbs3, Tsp3

Immunogen:

Recombinant Mouse Thrombospondin-3 protein (24-350AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A