

# Thbs3 Antibody



PACO30382

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse, Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000

**Protein Background:**

**Gene ID:**

Thbs3

**Uniprot**

Q05895

**Synonyms:**

Thrombospondin-3, Thbs3, Tsp3

**Immunogen:**


Recombinant Mouse Thrombospondin-3 protein (24-350AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



### Western blot

All lanes: Thbs3 antibody at 16 $\mu$ g/ml + HepG2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 105 kDa

Observed band size: 105 kDa

.