

# TK1 Antibody, FITC conjugated



PACO30418

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TK1

**Uniprot**

P04183

**Synonyms:**

Thymidine kinase, cytosolic (EC 2.7.1.21), TK1

**Immunogen:**

Recombinant Human Thymidine kinase, cytosolic protein (2-234AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A