

TMEM120B Antibody



PACO30466

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Necessary for efficient adipogenesis.

Gene ID:

TMEM120B

Uniprot

A0PK00

Synonyms:

Transmembrane protein 120B, TMEM120B


Immunogen:

Recombinant Human Transmembrane protein 120B protein (1-101AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO30466 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO30466 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO30466 at dilution of 1:100.