

TMEM154 Antibody



PACO30494

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

TMEM154

Uniprot

Q6P9G4

Synonyms:

Transmembrane protein 154, TMEM154


Immunogen:

Recombinant Human Transmembrane protein 154 protein (97-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: TMEM154 antibody at 6 μ g/ml + MCF-7 whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution


Predicted band size: 21 kDa

Observed band size: 21 kDa

.



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO30494 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer tissue using PACO30494 at dilution of 1:100.