

TMEM9 Antibody



PACO30654

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in intracellular transport.

Gene ID:

TMEM9

Uniprot

Q9P0T7

Synonyms:

Transmembrane protein 9 (Dermal papilla-derived protein 4), TMEM9, DERP4 TMEM9A

Immunogen:

Recombinant Human Transmembrane protein 9 protein (111-183AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO30654 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO30654 at dilution of 1:100.